STERILE	EO
---------	----

## **PBMV Balloon Catheter Set**

**(Percutaneous Balloon Mitral Valvuloplasty: PBMV)**

### **Instructions for Use**

**Doc. No. : SYM/TF-05-1902-03**

**Rev: V1.1**

**Revision Date: 2023.11.02**

# **CE2797**

Instructions for Use	( page 1~19)
Instruções de utilização	( página 20~39 )
Mode d'emploi	( page 40~59 )
Istruzioni per l'Uso	( pagina 60~79 )
Instrucciones de uso	( pagina 80~99 )
Handbuch für das Produkt	( seite 100~119 )
使用说明书	( page 120~136 )

# Contents

1. Device Description .....	1
2. Device Specifications and Size Parameters .....	2
3. Selection of PBMV Balloon Catheter .....	2
4. Intended Use .....	3
5. Indications .....	3
6. Contraindications .....	4
7. Potential Adverse Events .....	4
8. Precautions and Warnings .....	6
9. Directions for Use .....	7
10. Shelf Life .....	17
11. Stock, Supply and Disposal .....	17
12. Package and Contents .....	17
13. Disclaimer of Warranty and Limitation of Remedy .....	17

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
# 1. Device Description

The PBMV Balloon Catheter Set consists of PBMV Balloon Catheter and accessories such as stretching tube, dilator, guide wire, style, syringe and ruler.


- **Description of PBMV Balloon Catheter**

The PBMV balloon catheter is an instrument for dilating the stenotic mitral orifice in patients with mitral stenosis. This catheter consists of an inner tube, an outer tube, and a rubber balloon that were connected together by a four-way stopcock or a three-way stopcock. The elongated balloon catheter by the stretching tube is inserted into the left atrium over a specially designed guide wire. Then the distal part of the balloon is inflated by injecting 5-7mL diluted contrast medium from the injecting port and the guide wire is removed. Direct the balloon catheter to the mitral orifice by manipulating the stylet and extract the stylet 4-5cm will advance the balloon catheter into the left ventricle. With the distal portion of the balloon inflated, the catheter is pulled back to position the balloon at the mitral orifice, and the remaining diluted contrast medium is rapidly injected to inflate the whole balloon. As soon as the predetermined amount of diluted contrast medium is totally injected into the balloon, the balloon should be deflated immediately, which ends one dilation of the mitral valve. (See Fig. 1-2)

- 1. Shrunk waist    2. Balloon    3. Expanded orifice    4. Outer Tube
- 5. Four-way Stopcock    6. Vent Port    7. Connecting Tube
- 8. One-way Stopcock    9. Inner Tube    10. Cap

Fig. 1 The Structure of Balloon Catheter (four-way stopcock)

- 1. Shrunk waist    2. Balloon    3. Expanded orifice    4. Outer Tube
- 5. Three-way Stopcock    6. Cap    7. Inner Tube

Fig. 2 The Structure of Balloon Catheter (three-way stopcock)

- **Description of Accessories**

- Stretching tube: Stretching balloon catheter. It is made from stainless steel tube with O. D. of 1.2mm, and Length of 800mm.
- Dilator: A thick tube with pyramidal shape in the distal side and used to dilate the insertion area. It is made from Polyethylene(PE)with O. D. of 14 Fr, and Length of 700mm.
- Guide wire: Made from stainless steel with pyramidal shape in the distal part and is welded with a spring. It is used to guide the catheter into the mitral valve. O. D. : 0.025"; length: 1800mm.
- Stylet: Directing balloon to mitral valve. Made from stainless steel with "J" shape and is welded with a spring. O. D. : 0.038"; length: 800mm.
- Syringe: Used to inflate balloon. Specifications: 30mL.
- Ruler: Used to measure the dimensions of the Balloon; the material is acrylonitrile butadiene styrene(ABS). Specification: 50mm.

## 2. Device Specifications and Size Parameters

The specifications and size parameters of PBMV Balloon Catheter Set are shown in Table 1-2

Table 1 The Specifications of PBMV Balloon Catheter

Specification		Range of Inflation Diameter
Four-way stopcock	Three-way stopcock	
PBMV-30S	PBMV-30YS	28mm~30mm
PBMV-28S	PBMV-28YS	26mm~28mm
PBMV-26S	PBMV-26YS	24mm~26mm
PBMV-24S	PBMV-24YS	22mm~24mm
PBMV-22S	PBMV-22YS	20mm~22mm
PBMV-20S	PBMV-20YS	18mm~20mm

Table 2 The size of balloon catheter and auxiliary instruments

Name of Item	Size		
	Outer diameter	Length	Volume
Balloon Catheter	4mm(12Fr. )	70cm	_____
Stretching Tube	1.2mm	80cm	_____
Dilator	4.67mm(14Fr. )	70cm	_____
Guide wire	0.64mm(0.025")	180cm	_____
Stylet	0.96mm(0.038")	80cm	_____
Syringe	_____	_____	30mL
Ruler	_____	50mm	_____

## 3. Selection of PBMV Balloon Catheter

Please refer to the following instructions before selecting the size of the balloon catheter:

- A. Select the balloon diameter according to the Dai's formula:  
Reference size(mm)=height(cm)/10+10
- B. Use the reference size as the balloon diameter for the first inflation in patients with optimal indications.
- C. Use the reference size minus 2mm as the balloon diameter for the first inflation in patients with suboptimal indications.
- D. The three-way stopcock and four-way stopcock PBMV balloon catheter only different in the way of air expelled. No matter which PBMV balloon catheter, the performance of the same size balloon catheter is equivalent. The choice will be based on doctor's requirements. Doctors, accustomed to the use of PTCA balloons, are recommended to choose three-way stopcock PBMV balloon catheter. The others are recommended to choose four-way stopcock PBMV balloon catheter.

#### **4. Intended Use**

The PBMV Balloon Catheter Set is intended to dilate the mitral stenosis using percutaneous balloon mitral valvuloplasty(PBMV)technique, which offers minimal invasive therapy for patients with mitral stenosis.

#### **5. Indications**

- Symptomatic isolated mitral stenosis(moderate to severe), NYHA function classII~III.
- No deformity of the mitral valve, no severe valve calcification, no subvalvular abnormality and no marked valve rigidity detected by echocardiography. Mitral valve area less than or equal to 1. 5cm<sup>2</sup>, and no left atrial thrombus.
- Sinus rhythm demonstrated by electrocardiography.
- Mean left atrial pressure greater than 1. 46KPa or 11mmHg, mean mitral valve pressure gradient during diastole greater than or equal to 1. 06KPa or 8mmHg.
- As for patients with restenosis after surgical commissurotomy, atrial fibrillation, mitral valve calcification, mild mitral or aortic regurgitation, PBMV can also be performed but with caution.
- Patients with severe pulmonary hypertension, left-or right-sided heart failure, or with significant renal or hepatic disease.
- Patients unsuitable to accept valve replacement because of other conditions such as pregnancy or have contraindications for long-term anticoagulation.

##### **Note: Assessment of the Efficacy of PBMV**

**A. Mitral diastolic murmur disappears or markedly diminished, and NYHA function class elevate at least one grade.**

**B. PBMV is considered successful when mean left atrial pressure less than or equal to 1.46KPa or 11mmHg, mean mitral valve pressure gradient less than or equal to 1.06KPa or 8mmHg, optimal when mean mitral valve**

pressure gradient less than or equal to 6mmHg. Usually, cardiac output will increase and pulmonary resistance decrease.

**C. The procedure is considered optimal when mitral valve area more than or equal to 2cm<sup>2</sup> demonstrated by two-dimensional echocardiography. For patients with suboptimal indications, the procedure is considered successful when mitral valve area more than or equal to 1.5 cm<sup>2</sup>.**

**D. No serious complications.**

## **6. Contraindications**

- Patients with active rheumatic fever, history of systemic embolization, and severe dis-arrhythmia.
- Patients with severe deformity of mitral valve or subvalve structure.
- More than moderate mitral or aortic regurgitation.
- Patient with contraindications for atrial septal puncture, such as giant right atrium, aneurysmatic dilatation of the aortic root, significant displacement or heterotaxia of the heart, history of systemic embolization.

## **7. Potential Adverse Events**

### **1) Possible complications**

Possible complications associated with use of the PBMV balloon catheter include, but are not limited to

- Death
- Increase in valvular regurgitation.
- Vasovagal reaction.
- Cardiac dysrhythmias such as atrial fibrillation, heart block sinus bradycardia, and ventricular achy cardiac or fibrillation.
- Congestive heart failure or pulmonary edema.
- Vascular with or without the need for surgical intervention.
- Injury to, or perforation of, the myocardial wall with or without cardiac tamponade.
- Damage to the mitral valve or ruptured chordae tendinitis.
- Thrombosis; Embolism cause by thrombus, valve fragment, calcification apparatus, or air bubbles.
- Myocardial infarction.
- Transient ischemic attack; transient visual defect.
- Blood loss requiring fluid replacement or transfusion of blood components; Hematoma.
- Atrial septal defect, transient or persistent, with or without hemodynamic compromise.

- Infection.
- Allergy or anaphylactic reaction to dilute contrast medium or device components.
- Decreased blood pressure when blood flow across the valve orifice is temporarily obstructed by the inflated balloon.
- Respiratory arrest.
- Renal failure.

## 2) Failure of the balloon to deflate

If the balloon is subjected to twisting forces beyond those which it is designed to withstand, a small opening in the inner balloon layer may develop, thus allowing leakage of the radiopaque dilute contrast medium into the space between the inner and outer latex layers of the balloon. At this point, the outer wall of the balloon is the only layer which will expand (see Fig3).


Fig. 3 The Structure of Balloon

Should this occur, the clinician would observe that deflation of the balloon following valvuloplasty becomes impossible. This is due to the continued flow of dilute contrast medium from the center of the balloon into the space between the two layers. The opening between the layers acts as a one-way valve and obstructs the return flow of dilute contrast medium as the inner layer contracts, thus preventing the dilute contrast medium from being withdrawn when negative pressure is exerted on the syringe.

If the balloon cannot be deflated, it must immediately be withdrawn from the valve orifice into the left atrium then carefully pulled just to the atrial septal wall(see Fig. 4).


Fig. 4 Withdraw the balloon

The unique design of the PBMV balloon catheter features two small holes (safety hole) through which the dilute contrast medium will gradually escape.

They are located 180° apart in the outer latex layer of the waist of the balloon. After withdrawing the balloon to the atrial septal wall, the user should then extend the balloon by inserting the guide wire and balloon stretching tube. The guide wire must extend beyond the catheter tip (see figure5).


Fig. 5 Safety hole

This action will facilitate and accelerate the flow of the dilute contrast medium from the balloon, which will take approximately 10 minutes. If deflation is not complete within this time, wait several more minutes until the balloon has deflated spontaneously further.

At this time, ensure that the heparin has been administered in adequate amounts to prevent thrombus formation; an additional dose may be required. The catheter may be withdrawn across the atrial septum when the balloon diameter is 15mm or less, as the balloon will be soft and pliable due to the outer layer has continued to expand.

**Caution: In this situation, do not use normal saline or other fluid to inflate the balloon. Do not attempt to expand the balloon to the point of rupture, as several times the maximum inflation volume is required for rupture to occur.**

## 8. Precautions and Warnings

- Read the instruction carefully before using the products, make sure if there is damage of the balloon catheter and its accessories. Don't use the product when the package is damaged. Please use the product before expiration date.
- Do not reuse. Discard after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. All parts are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused. Accordingly, SYM will not be responsible for any direct, incidental or consequential damages resulting from reuse of the product.
- The products should be used only by doctors who are trained and qualified in PBMV techniques. Doctors who use this product should keep the indications and contraindications in mind, and select the balloon diameter according to 3.
- During the procedure, if the balloon can't be deflated, the operators should make

his best to withdraw the balloon into left atrium. The cause of this situation is rupture of the inner layer of the balloon resulting from manipulation mistake. To avoid this, never push the inner tube during balloon preparation unless a stretching tube and guide wire has been placed in it. When the balloon is withdrawn into the left atrium, stretch the balloon to force the liquid out of the balloon.

- During the procedure, patients should be heparinized (100 u/Kg) and all the instruments be flushed with heparinized saline. Never touch the balloon segment of the catheter with hand.
- If the valve is severely sclerotic and/or highly calcified or if severe subvalvular stenosis is present, balloon malfunction or tearing of the mesh between the outer and inner balloon membranes could occur. If at the time of inflating the balloon, an unusual shape is observed, immediately deflate the balloon, withdraw it and replace it with a new one.
- This PBMV procedure must be performed in facilities where cardiac surgeons and their surgical back-up are available in-house.
- This product contains nature rubber latex, which may cause such allergic reactions as: itching, redness, urticaria, swelling, fever, dyspnea, asthma-like symptoms, hypotension, shock, etc. Should one or more of these symptoms be observed in a patient, stop the procedure immediately and implement appropriate treatment.
- Examine the products with care, prior to use. Always keep extra sets of the products so that replacement can be made if necessary.
- Use in Specific population

Pregnancy - Care must be taken to minimize the radiation exposure to the fetus and the mother.

Nursing mothers - There has been no quantitative assessment of the presence of leachable in breast milk.

## 9. Directions for Use

The recommended use of this product is shown as following:

- **Pre-procedure**

- 1) Flush the inner tube with heparinized saline (see Fig. 6).**


Fig. 6 Flush the inner tube

## 2) Venting the balloon Catheter

- a) **Four-way stopcock type** Open the injecting port, keep the vent port close. Inject saline into the catheter slowly while the balloon is placed at a lower position. Suck the saline out while the balloon is Inflated. If necessary repeat this procedure until all air is removed. Open the vent port and inject saline till the saline flow out from the vent port. Close the injecting port and the vent port (see fig. 7).


Fig. 7 Venting the four-way stopcock balloon Catheter

- b) **Three-way stopcock type** Open the injecting port. Inject saline into the catheter slowly while the balloon is placed at a lower position. Suck the saline out while the balloon is Inflated. If necessary repeat this procedure until all air is removed. Then close the injecting port (see fig. 8).


Fig. 8 Venting the three-way stopcock balloon Catheter

**Caution: Unless air is expelled from the balloon lumen, the balloon may fail to inflate up to its specified diameter. If the balloon ruptures, air embolization may result. To avoid this discharge the air from the balloon lumen completely.**

## 3) Pretesting the balloon

Fill the syringe with the amount of dilute contrast medium corresponding to the maximal balloon diameter. Inflate the balloon and measure the balloon diameter with the ruler. Correct the difference between the actual diameter and the

designed diameter by increasing or decreasing the injected amount of fluid.

It is recommended to check the sequence of balloon inflation and the structure of the balloon catheter whether or not fulfilling the designing standards. The normal sequence of balloon inflation is: first, the distal part, then the proximal part, and finally the waist (see fig. 9).


Fig. 9 Pretesting the balloon


Fig. 10 Flushing the Stretching Tube

#### 4) Flushing the Stretching Tube

Flush the stretching tube with heparinized saline (see fig. 10).

#### 5) Stretching the Balloon Catheter

Insert the stretching tube into the inner tube and lock it with the inner tube. Push both the inner tube and the stretching tube until the pin of the inner tube is locked in the slot. By doing so, the balloon segment of the catheter is stretched (see fig. 11 and fig. 12).


Fig. 11 Stretching the Balloon Catheter 1



Fig. 12 Stretching the Balloon Catheter 2

**Caution:** After the balloon has been elongated, do not disconnect the metal hubs from each other or the four-way stopcock unless the guide wire is in place and fully protrudes from the catheter tip. Do not advance the catheter inner tube into the four-way stopcock unless the stretching tube has been inserted into the catheter inner tube and Luer lock at the hubs. Should this caution be disregarded, the catheter could bend at an acute angle resulting in damage to the device or difficulty in carrying out the procedure. Never use a kinked balloon.

#### 6) Preparing the Stylet

The stylet is made of spring wires. Its shape can be pre-shaped by operator so as to facilitate the balloon catheter advancing into the left ventricle, according to the size of the left atrium(see fig. 13).

Used for moderate sized left atrium (35~45mm)

Used for larger left atrium (>50mm)

Used for smaller left atrium( $\leq 35$ mm)

Fig. 13 Preparing the Stylet

**Caution:**

a. During manipulation of the stylet (to direct the balloon), it is very important to hold the four-way stopcock rather than the inner tube (proximal end of the catheter).

If the inner tube is held, the lumen of the catheter may become twisted, thus occluding the lumen and blocking subsequent advancement or withdrawal of the guide wire or stylet. Do not twist the catheter lumen in this manner as it may adversely affect balloon performance.

b. A spring coil is wound around a solid core in the stylet. The stylet may be manually shaped, but repeated forced bending of the stylet may cause damage to the core wire resulting in separation of the wire, and the spring wire may uncoil. Shape the stylet, taking care not to kink or separate the core wire.

**7) Puncturing the interatrial Septum**

Puncture the interatrial septum by a standard method.

- **Operation**

**1) Inserting the Balloon Catheter Into the Left Atrium**

**a. Dilating the Puncture Site of Atrial Septum**

Insert the dilator over the guide wire to dilate the groin entry site, the femoral vein and the atrial septal

puncture site. After doing so, remove the dilator (see fig. 14).


Fig. 14 Dilating the Puncture Site of Atrial Septum

**b. Inserting Balloon Catheter**

To facilitate passing through the skin and soft tissue at groin entry site, the balloon segment can be covered with contrast medium or silicon oil. Insert the elongated balloon over the guide wire.

Advance the balloon catheter over the guide wire until a slight resistance is encountered at the atrial septum. At this site, the manipulation should be gentle. If the passage is not easy, adjust the catheter tip will enable the balloon catheter pass into the left atrium.

The stretching tube and inner tube must be locked with the four-way stopcock before the balloon catheter passes through the atrial septum.

It is important to follow the above steps when inserting the balloon catheter into the left atrium (see fig. 15).


Fig. 15 Inserting Balloon Catheter into left atrium

**c. Extracting Stretching Tube**

When the greatest portion of the balloon has entered into the left atrium and the tip of the balloon nears the roof of the left atrium, pull the stretching tube approximately 2-3cm. This will make the distal portion of the balloon more elastic (see fig. 16).


Fig. 16 Extracting Stretching Tube

**d. Advancing Balloon Catheter into the Left Atrium**

After pulling the stretching tube 2-3cm, push the balloon catheter over the guide wire about 5cm, thus passing the balloon segment into the left atrium completely (see fig. 17).


Fig. 17 Advancing Balloon Catheter into the Left Atrium

### e. Extracting the Inner Tube

Unlock the inner tube from the four-way stopcock carefully and pull the inner tube out until resistance is felt.

**Caution: Be careful to unlock the inner tube, and not to pull the balloon catheter out of the left atrium when extracting the inner tube (see fig. 18).**


Fig. 18 Extracting the Inner Tube

### f. Advancing the Balloon Catheter into the Left Atrium

After the balloon segment has been placed in the left atrium, advance the balloon catheter further over the guide wire so that there will be sufficient length of the catheter placed in the left atrium. This also makes the catheter have a loop that will enable the catheter cross the mitral orifice easily by manipulating the stylet.

**Caution: During the above steps, it is important not to pull the guide wire out of the left atrium (see fig. 19).**


Fig. 19 Advancing the Balloon Catheter into the Left Atrium

## 2) Crossing the Mitral Orifice

### a. The Counterclockwise Rotation Method:

Inflate the distal portion of the balloon partially with a small amount of diluted contrast medium, then remove the guide wire and stretching tube together. Insert stylet into the inner tube and rotate the stylet counterclockwise. This will direct the balloon tip toward the mitral orifice; When the balloon tip is near the mitral orifice the catheter tip will have a “wood-pecking” movement. At this time, pulling the stylet

out will make the balloon move forward and enter into the left ventricle. If the mitral orifice is too small, the distal portion of the balloon may not be inflated to aid the balloon crossing the mitral orifice (see fig. 20).


Fig. 20 The Counterclockwise Rotation Method

**b. The Clockwise Rotation Method**

This method is useful when the left atrium is large (>55mm). Advance the balloon catheter 5~6cm into the left atrium (see fig. 21).


Fig. 21 The Clockwise Rotation Method 1

Rotate the stylet clockwise and make the catheter a large loop. This will bring the balloon tip toward the mitral orifice (see fig. 22).


Fig. 22 The Clockwise Rotation Method 2

With the stylet held firm in place, the balloon catheter is advanced, allowing the balloon move into the left ventricle along the inferior of the left atrium (see fig. 23).


Fig. 23 The Clockwise Rotation Method 3

### 3)Inflating the Balloon

a. Confirm the amount of diluted contrast medium in the syringe can inflate the balloon to the predetermined diameter. Close the vent port. With the distal portion of the balloon partially inflated balloon catheter gently to fix the balloon at the mitral orifice.

**Caution: Do not pull the catheter too strong or too light, otherwise it will disposition the balloon and affect the efficacy (see fig. 24).**


Fig. 24 Inflating the Balloon

b. When the balloon is positioned at the mitral orifice, inject the remaining diluted contrast medium rapidly to inflate the whole balloon and dilate the valve (a). As soon as the balloon is fully expanded, the balloon should be deflated and be placed in the atrium (b) (see Fig. 25).

#### **Caution:**

**a. The duration of the balloon inflation should be kept within 5 seconds so that**

mitral orifice obstruction by the inflated balloon will not result in severe complications.

b. After inflation, pressures in left atrium and ventricle should be measured, and echocardiography be performed. If the hemodynamics and the mitral valve opening are satisfactory, the balloon catheter can be removed and the procedure discontinued.


Fig. 25 Balloon be deflated and placed in the atrium

#### 4) Extracting the Balloon Catheter

First, insert the guide wire into the stretching tube, then insert the guide wire and the stretching tube together into the inner tube. When the stretching tube and the guide wire near the tip of the balloon, push the guide wire out of the catheter into the left atrium ( see Fig. 12). Lock the stretching tube with the inner tube together. Under fluoroscopic guidance, push the inner tube and the stretching tube together all the way and lock them in the slot. Then pull the stretched balloon catheter and accessories as a unit into right atrium or inferior vena cava.

##### Caution:

a. It is important to stretch the balloon catheter under fluoroscopic guidance, and make sure there is sufficient length of guide wire in the left atrium so as not to perforate the left atrium.

b. The soft tip of the guide wire must protrude from the catheter tip all the time during the above steps.

## **10. Shelf Life**

The PBMV Balloon Catheter Set is sterilized by ethylene oxide. The shelf life is 24 months. Please use the device before the “Use By” date printed on the label.

## **11. Stock, Supply and Disposal**

- The balloon catheter and its accessories are disposable. Please contact our distributor or manufacturer if any abnormality of this product is detected and do not use it.
- Do not use the product beyond the expiration date.
- In order to provide maximum protection of the PBMV balloon catheter, it must be stored in a dry and cool place free from contamination. It is recommended the stored catheters be checked regularly in order to ensure the normal operation and to avoid exceeding the “Use By” date.

## **12. Package and Contents**

The product package consists of an outer paper box, an inner paper box and the sterilized pack. The package is marked by product model, sterilization date, expiration date, serial number, sterilization indicator.

## **13. Disclaimer of Warranty and Limitation of Remedy**

- SYM will replace any defective the PBMV Balloon Catheter and set contents (e. g. , damaged product), free of charge. In the event of a product complaint, the user is requested to return the PBPV Balloon Catheter and packaging materials back to the distributor or manufacturer so that the cause of the complaint may be fully investigated.
- SYM is not responsible for any damage caused due to improper handling, operation and storage, including use after the expiry date printed on the label.
- SYM is not responsible for any damage to the PBMV Balloon Catheter or its accessories caused by transportation, handling and storage at the health care facility, whether physical damage or human injury.
- SYM is not responsible for any damage caused by reprocessing or reuse of the PBMV Balloon Catheter or its accessories, whether physical damage or human injury.
- SYM is not responsible for any impairment to structural integrity and/or function which may arise from the use of the PBMV Balloon Catheter or its accessories with other manufacturer’s products. Do not use the PBMV Balloon Catheter or its accessories with other manufacturer’s products.
- SYM is only responsible for any damage to patients and users caused by the intended use of PBMV balloon catheter set, SYM won’t responsible for any damage to patients and users caused by unexpended use of the PBMV balloon catheter set.

SYMBOL	SYMBOL DEFINITIONS
	Manufacturer
	Authorized representative in the European Community/ European Union
	Country of manufacture
	Model number
	Use-by date
	Batch code
	Serial number
	Sterilized using ethylene oxide
	Do not re-sterilize
	Do not use if package is damaged and consult instructions for use
	Do not reuse
	Consult instructions for use or consult electronic instructions for use
	Caution
	Contains or presence of natural rubber latex
	Medical device
	Unique device identifier
	Single sterile barrier system

	Single sterile barrier system with protective packaging outside
	Fragile, handle with care
	Keep away from sunlight
	Keep dry
<b>CONTENTS</b>	Contents of the package
<b>LENGTH</b>	Usable length of the device
<b>CE2797</b>	CE marking 2797: Identification number of notified body

**Manufacturer:**

**Shenzhen Shineyard Medical Device Company Limited**

3/F, Changfeng Industrial Block, No.3 Liuxian Road, Xin'an, Bao'an District, Shenzhen, Guangdong, 518000, China

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SUNGO Europe B.V.